DOPPLER EFFECT & STELLAR SPECTRA

8B20.20

Stellar Spectra







Equipment:

- Balls with Dowel Arms (2)

Concept:

The white ball represents the center of a galaxy. The colored balls 1. Powered Rotator represent stars revolving around the galactic center. The amount of 2. Large Styrofoam Ball red/blue seen represents the degree to which the spectrum is Doppler 3. Small Red/Blue Styrofoam shifted. A star's spectrum is shifted toward the blue as the star approaches the observer because its wavelength is compressed. As the star recedes

from the observer, the spectrum shifts toward the *red* since its wavelength is extended. A useful mnemonic device is "red for recede".

The Doppler shifted frequency for light with observer and source moving *toward* and (*away from*) each other with relative speed v are given, respectively by

$$f_{\text{observer}} = \sqrt{\frac{1 + \frac{v}{c}}{1 - \frac{v}{c}}} f_{\text{source}} \text{ and } f_{\text{observer}} = \sqrt{\frac{1 - \frac{v}{c}}{1 + \frac{v}{c}}} f_{\text{source}}$$

The effect clearly becomes more pronounced as the relative speed v, approaches the speed of light, c.

Procedure:

- 1. Verify that the white ball is securely attached to the rotator base and the colored balls are arranged as shown in the top-left picture.
- 2. Toggle the rotation direction switch to counter-clockwise (CCW).
- 3. Turn the rotation speed knob clockwise to turn the powered rotator on and adjust its speed.
- 4. Use the switch on the power strip (not shown) to turn the motor on and of f to avoid being hit by the revolving balls.