OHM'S LAW DC Circuits Electricity & Magnetism

Concept:

This demonstration verifies the relation V = IR, where V is potential 1. DMM - Voltmeter difference, I is current, and R the resistance. Some texts quibble with this 2. Battery Eliminator (1.5, 3, 4.5, 4.5)as simply a definition of resistance, and claim that the correct statement of Ohm's Law is $J = \sigma E$, where J is the current density, σ is conductivity, and **E** the electric field. Perhaps more importantly, only certain materials, called ohmic, obey these relations. Examples of nonohmic materials abound in modern microelectronics; many of the semiconductor components of a computer are prime examples.



5F10.10



Equipment:

- 6, 9, 12 V)
- 3. DMM Ammeter
- 4. (3) 3 ft, Red Banana-Banana Cables
- Variable Resistor
- 6. (2) 3ft, Black Banana-Banana Cables
- 7. (2) Alligator Clips

Procedure:

- 1. Using one of the DMMs as an ohmmeter, verify that the variable resistor is set to 10 Ω and the circuit is wired as shown in the picture above.
- 2. Tell the class that the resistor is set to ~ 10 Ω and the voltage is set to ~ 1.5 V.
- 3. Ask the class to use ohm's law to predict what the current should be.
- 4. Turn the battery eliminator on and reveal the measured current and voltage.
- 5. Repeat steps 2-4 after varying the voltage.
- 6. Repeat steps 2-4 after varying the resistance (be careful not to exceed 1 A as this could blow the battery eliminator's fuse).