JUMPING RING

5K20.30

Launch Force

(**F**)

11:

Eddy Currents

→ B
_

Electricity and Magnetism

Electromagnetic Induction



Concept:

This is a crowd pleaser with subtle physics. The upward force on the ring • requires that the coil's B-field have a *radial* component, and this arises only because the field fringes as shown in the diagram above. The upward launching force on the ring is given by:

$$\mathbf{F} = I\mathbf{L} \times \mathbf{B}$$

where *I* is the induced current in the ring with direction given by Lenz's Law, **L** is the current segment length in the direction of the current, and \mathbf{B}_{r} is the radial component of the coil's magnetic field. Note that there is also an *axial* component of **B** (not shown), which exerts a *compressive* force on



Magnetic Field

(B)

Current

4-----

Equipment:

- Ring Launcher Apparatus
- (2) Large Aluminum Rings
- Induction Coil and Bulb
- Copper Ring
- Split Aluminum Ring
- Small Aluminum Ring
- Liquid Nitrogen*, tongs and Dewar (upon request)

the ring but with unseen effect because of the ring's stiffness. The reason this demonstration works with AC supplied to the coil is due to a phase difference between the supplied AC and the current induced in the ring (see Am. J. Phys. **79**, 353 (2011) and Am. J. Phys. **68**, 238 (2000)).

Procedure:

- 1. Verify that the ring launcher is plugged in and the green light is on.
- 2. Slide one of the rings onto the center rod.
- 3. Press the "launch" switch to activate the ring launcher. Don't depress the switch for more than 1 sec!
- * Be aware of low ceilings. Liquid nitrogen cooled rings may launch substantially higher than 1 meter.

	Large Aluminum Ring	Small Aluminum Ring	Copper Ring	Split Aluminum Ring	Induction Coil & Bulb
Action	Launches ~ 2 m	Launches ~ 1 m	Launches ~ 1 m	None	Lights up